

# Recombinant Aspergillus niger Aspergillopepsin-2

Catalog No: #AP74083

Package Size: #AP74083-1 20ug #AP74083-2 100ug #AP74083-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Recombinant Aspergillus niger Aspergillopepsin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:60-98aaSequence Info:Partial
Other Names	Acid protease A Aspergillopepsin II Proctase A
Accession No.	P24665
Uniprot	P24665
Calculated MW	19.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	EEYSSNWAGAVLIGDGYTKVTGEFTVPSVSAGSSGSSGY
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## References

"The gene and deduced protein sequences of the zymogen of Aspergillus niger acid proteinase A."Inoue H., Kimura T., Makabe O., Takahashi K.J.  
Biol. Chem. 266:19484-19489(1991)Research Topic:Microbiology

---

Note: This product is for in vitro research use only