

VPS53 Antibody

Catalog No: #24678




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	VPS53 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB ICC
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide
Immunogen Description	Raised against a 18 amino acid peptide from near the carboxy terminus of human VPS53.
Target Name	VPS53
Other Names	Vacuolar protein sorting 53, HCCS1
Accession No.	Swiss-Prot:Q5VIR6Gene ID:55275
Uniprot	Q5VIR6
GeneID	55275;
Concentration	1mg/ml
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of VPS53 in 293 cell lysate with VPS53 antibody at (A) 0.5 and (B) 1 ug/mL.



Immunocytochemistry of VPS53 in 293 cells with VPS53 antibody at 5 ug/mL.

Background

The sorting of acid hydrolases to lysosomes rely on mannose 6-phosphate receptors that cycle between the trans-Golgi network (TGN) and endosomes. The maintenance of this cycle requires the function of the mammalian Golgi-associated retrograde protein (GARP) complex which is composed of three subunits: VPS52, VPS53, and VPS54. Depletion of any of these three proteins, such as by RNAi, impairs the retrograde transport of multiple TGN proteins. VPS53 was identified as an HIV dependency factor (HDF) and plays a role in viral entry to the cell, suggesting that VPS53 may be an important drug target in HIV treatment. At least five isoforms of VPS53 are known to exist.

Note: This product is for in vitro research use only