

MLANA MART1 Antibody FITC Conjugated

Catalog No: #C00036F



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | MLANA MART1 Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF |
| Species Reactivity | Hu |
| Immunogen Description | KLH conjugated synthetic peptide aa 50-100 118 derived from human MLANA MART1 |
| Conjugates | FITC |
| Target Name | MLANA MART1 |
| Other Names | MART1; MART-1; Melanoma antigen recognized by T-cells 1; Antigen LB39-AA; Antigen SK29-AA; Protein Melan-A; MLANA |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Intracellular |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Background

Involved in melanosome biogenesis by ensuring the stability of GPR143. Plays a vital role in the expression, stability, trafficking, and processing of melanocyte protein PMEL, which is critical to the formation of stage II melanosomes.

Note: This product is for in vitro research use only