

beta Actin Antibody HRP Conjugated

Catalog No: #C00044H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | beta Actin Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | WB IHC-P IHC-F |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | MAP conjugated synthetic peptide aa 1-50 375 derived from human beta Actin |
| Conjugates | HRP |
| Target Name | beta-Actin |
| Other Names | BRWS1; PS1TP5BP1; Actin, cytoplasmic 1; Beta-actin; ACTB |
| Accession No. | Swiss-Prot#P60709NCBI Gene ID60 |
| Uniprot | P60709 |
| GeneID | 60; |
| Excitation Emission | N A |
| Cell Localization | Cytoplasm |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:10000-1:50000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

Note: This product is for in vitro research use only