

NFkB Inducing Kinase NIK Antibody HRP Conjugated

Catalog No: #C00056H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NFkB Inducing Kinase NIK Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 900-947 947 derived from human NIK
Conjugates	HRP
Target Name	NFkB Inducing Kinase NIK
Other Names	HS; NIK; HSNIK; FTDCR1B; Mitogen-activated protein kinase kinase kinase 14; NF-kappa-beta-inducing kinase; Serine threonine-protein kinase NIK; MAP3K14
Accession No.	Swiss-Prot#Q99558NCBI Gene ID9020
Uniprot	Q99558
GeneID	9020;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Lymphotoxin beta-activated kinase which seems to be exclusively involved in the activation of NF-kappa-B and its transcriptional activity. Promotes proteolytic processing of NFkB2 P1, which leads to activation of NF-kappa-B via the non-canonical pathway. Could act in a receptor-selective manner.

Note: This product is for in vitro research use only