

Hepatitis E Virus ORF3 Antibody HRP Conjugated

Catalog No: #C00152H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatitis E Virus ORF3 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	HEV
Crossing Reactivity	HEV
Immunogen Description	KLH conjugated synthetic peptide aa 85-104 114 derived from HEV ORF3
Conjugates	HRP
Target Name	Hepatitis E Virus ORF3
Other Names	structural protein; immunogenic protein; ORF-3; ORF 3; Protein ORF3; pORF3; phosphoprotein; ORF3 [Hepatitis E virus]; ORF3_HEVHY.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Hepatitis E Virus ORF3

Note: This product is for in vitro research use only