

bFGF Antibody HRP Conjugated

Catalog No: #C00157H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	bFGF Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 20-60 155 derived from mouse bFGF
Conjugates	HRP
Target Name	bFGF
Other Names	Fgfb; bFGF; Fgf-2; Fibroblast growth factor 2; Basic fibroblast growth factor; Heparin-binding growth factor 2; HBGF-2; Fgf2
Accession No.	Swiss-Prot#P15655NCBI Gene ID14173
Uniprot	P15655
GeneID	14173;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration. Functions as potent mitogen in vitro (By similarity).

Note: This product is for in vitro research use only