

## FGF1 Antibody HRP Conjugated

Catalog No: #C00166H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FGF1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from mouse FGF1
Conjugates	HRP
Target Name	FGF1
Other Names	Fam; Fgfa; Dffrx; Fgf-1; Fibroblast growth factor 1; Acidic fibroblast growth factor; aFGF; Heparin-binding growth factor 1; HBGF-1; Fgf1
Accession No.	Swiss-Prot#P61148NCBI Gene ID14164
Uniprot	P61148
GeneID	14164;
Excitation Emission	N A
Cell Localization	Cytoplasm, Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration. Functions as potent mitogen in vitro (By similarity).

Note: This product is for in vitro research use only