

BSA Antibody HRP Conjugated

Catalog No: #C00211H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BSA Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Cow
Immunogen Description	BSA
Conjugates	HRP
Target Name	BSA
Other Names	Serum albumin; BSA; ALB
Accession No.	Swiss-Prot#P02769NCBI Gene ID280717
Uniprot	P02769
GeneID	280717;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 8% of all plasma zinc.

Note: This product is for in vitro research use only