

Mesothelin Antibody FITC Conjugated

Catalog No: #C00214F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mesothelin Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 130-180 630 derived from human Mesolin
Conjugates	FITC
Target Name	Mesothelin
Other Names	MPF; SMRP; Mesothelin; CAK1 antigen; Pre-pro-megakaryocyte-potentiating factor; MSLN
Accession No.	Swiss-Prot#Q13421NCBI Gene ID10232
Uniprot	Q13421
GeneID	10232;
Excitation Emission	494nm 518nm
Cell Localization	Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Background

Membrane-anchored forms may play a role in cellular adhesion. Megakaryocyte-potentiating factor (MPF) potentiates megakaryocyte colony formation in vitro.

Note: This product is for in vitro research use only