

General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit

Catalog No: #EK20168

Package Size: #EK20168-1 48T #EK20168-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit
Applications	ELISA
Species Reactivity	General

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.118ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only