

## FITC anti-human CD282 (TLR2)

Catalog No: #26656



Package Size: #26656-1 25 tests #26656-2 100 tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FITC anti-human CD282 (TLR2)
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1
Conjugates	FITC
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

10 µl/test

## Background

The monoclonal antibody reacts with mouse and human Toll-like receptor 2 (TLR2). Mouse TLR2 is expressed by the myeloid lineage, including macrophage and dendritic cells in splenocytes and the RAW264.7 cell line. To date, at least twelve members of

Note: This product is for in vitro research use only