

FITC anti-human CD33

Catalog No: #26660



Package Size: #26660-1 25 tests #26660-2 50 tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FITC anti-human CD33
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1
Conjugates	FITC
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

10 µl/test

Background

The antibody reacts with CD33 antigen, a 67 kDa type I transmembrane glycoprotein present on myeloid progenitors, monocytes and granulocytes. CD33 is absent on lymphocytes, platelets, erythrocytes, hematopoietic stem cells and non-hematopoietic system.. CD33 an

Note: This product is for in vitro research use only