

FITC anti-human CD45RB

Catalog No: #26666



Package Size: #26666-1 25 tests #26666-2 50 tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FITC anti-human CD45RB
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1, κ
Conjugates	FITC
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

10 μ l/test

Background

The monoclonal antibody 4APD7/26 recognizes CD45RB, which is an isoform of the Leukocyte Common Antigen (LCA) also known as CD45. CD45RB is a protein tyrosine phosphatase receptor that is present on most hematopoietic cells and absent on non-hematopoietic

Note: This product is for in vitro research use only