

FITC/PE/PerCP/APC anti-human CD3/CD8/ CD45/CD4

Catalog No: #26703

Package Size: #26703-1 20 tests #26703-2 50 tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FITC/PE/PerCP/APC anti-human CD3/CD8/ CD45/CD4
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG2a, κ/Mouse IgG1/Mouse IgG2a/Mouse IgG2b, κ
Conjugates	FITC/PE/PerCP/APC
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

20 µl/test

Background

FITC/PE/PerCP/APC anti-human CD3/CD8/ CD45/CD4

Note: This product is for in vitro research use only