

PE anti-human CD18

Catalog No: #26715



Package Size: #26715-1 25 tests #26715-2 50 tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-human CD18
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1
Conjugates	PE
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

10 µl/test

Background

The antibody reacts with CD18, the 95 kDa beta chain component of leukocyte function associated antigen-1 (LFA-1). CD18 is expressed by all peripheral blood leukocytes. CD18 is a leukocyte adhesion receptor that is essential for cell-to-cell contact in

Note: This product is for in vitro research use only