

PE anti-human CD21

Catalog No: #26719



Package Size: #26719-1 25 tests #26719-2 50 tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-human CD21
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1
Conjugates	PE
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

10 µl/test

Background

The antibody reacts with CD21 protein, a 145 kDa receptor for C3d complement fragment and for Epstein-Barr virus (EBV), expressed on all mature B lymphocytes from slg+ stage, then lost upon activation. It also present on follicular dendritic cells and some epit

Note: This product is for in vitro research use only