

PE anti-human CD279 (PD-1)

Catalog No: #26726



Package Size: #26726-1 25 tests #26726-2 100 tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-human CD279 (PD-1)
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1, κ
Conjugates	PE
Formulation	PBS, pH 7.2, with 0.1% gelatin, 0.2% BSA,0.09% sodium azide
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

5 µl/test

Background

The monoclonal antibody reacts with the human PD-1 (programmed death-1), a 55 kDa member of the immunoglobulin superfamily. PD-1 contains the immunoreceptor tyrosine-based inhibitory motif (ITIM) and plays a key role in peripheral tolerance and aut

Note: This product is for in vitro research use only