

## PE anti-human CD73

Catalog No: #26760



Package Size: #26760-1 25 tests #26760-2 100 tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE anti-human CD73
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1, $\kappa$
Conjugates	PE
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

5  $\mu$ l/test

## Background

CD73 is a 70 kD glycosylphosphatidylinositol (GPI)-linked 5'-nucleotidase, which is also known as ecto-5'-nucleotidase. It converts adenosine monophosphate (AMP) to adenosine. CD73 is expressed on subsets of T and B cells, mesenchymal stem cells, follicular dendritic cells, endothelial cells, and epithelial cells. It has been reported that CD73 costimulates T cell activation, and mediates adhesion of lymphocytes to follicular dendritic cells and endothelial cells.

Note: This product is for in vitro research use only