

PE anti-human CD9

Catalog No: #26766



Package Size: #26766-1 25 tests #26766-2 50 tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PE anti-human CD9
Host Species	Mouse
Clonality	Monoclonal
Applications	FC
Species Reactivity	Human
Immunogen Type	Mouse IgG1, κ
Conjugates	PE
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

10 µl/test

Background

The antibody reacts with CD9 (p24), a member of the tetraspan (TM4SF) family with 24 kDa MW, expressed on platelets and weakly on B-cells. It also expressed on eosinophils, basophils, endothelial and epithelial cells. CD9 antigen modulates cell adhesion, mig

Note: This product is for in vitro research use only