

Nestin Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90003PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nestin Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Ms Rat
Immunogen Type	681-750/1893
Conjugates	PE-Cy5.5
Target Name	Nestin
Other Names	Nestin
Accession No.	Uniprot:P21263Gene ID:25491
Uniprot	P21263
GeneID	25491
Calculated MW	209
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only