

# Amylin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90029Bio



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Amylin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-37/89
Conjugates	Bio
Target Name	Amylin
Other Names	Amylin
Accession No.	Uniprot:P10997Gene ID:3375
Uniprot	P10997
GeneID	3375
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only