

Neurokinin 1 Receptor Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90034PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurokinin 1 Receptor Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2-50/407
Conjugates	PE-Cy7
Target Name	Neurokinin 1 Receptor
Other Names	Neurokinin 1 Receptor
Accession No.	Uniprot:P25103Gene ID:6869
Uniprot	P25103
GeneID	6869
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only