

## CD19 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90043PE-Cy5.5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD19 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	485-556/556
Conjugates	PE-Cy5.5
Target Name	CD19
Other Names	CD19
Accession No.	Uniprot:P15391Gene ID:930
Uniprot	P15391
GeneID	930
Calculated MW	59
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only