

CD31 Polyclonal Antibody Cy5.5 Conjugated

Catalog No: #C90103Cy5.5



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD31 Polyclonal Antibody Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	681-738/738
Conjugates	Cy5.5
Target Name	CD31
Other Names	CD31
Accession No.	Uniprot:P16284Gene ID:5175
Uniprot	P16284
GeneID	5175
Calculated MW	78
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only