

GLP-2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90109Bio



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GLP-2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	146-165/180
Conjugates	Bio
Target Name	GLP-2
Other Names	GLP-2
Accession No.	Uniprot:P01275Gene ID:2641
Uniprot	P01275
GeneID	2641
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only