

E1 glycoprotein Polyclonal Antibody Bio Conjugated

Catalog No: #C90148Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E1 glycoprotein Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Rubellavirus
Conjugates	Bio
Target Name	E1 glycoprotein
Other Names	E1 glycoprotein
Accession No.	Uniprot:P08563Gene ID:
Uniprot	P08563
Calculated MW	52
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only