

MADH7/Smad7 Polyclonal Antibody FITC Conjugated

Catalog No: #C90237FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MADH7/Smad7 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/426
Conjugates	FITC
Target Name	MADH7/Smad7
Other Names	MADH7/Smad7
Accession No.	Uniprot:O15105Gene ID:4092
Uniprot	O15105
GeneID	4092
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only