

# Collagen I Polyclonal Antibody FITC Conjugated

Catalog No: #C90242FITC



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Collagen I Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	Collagen I
Other Names	Collagen I
Accession No.	Uniprot:P02452Gene ID:1277
Uniprot	P02452
GeneID	1277
Calculated MW	27/95kDa
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only