

CD3E Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90306PE-Cy7



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD3E Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Conjugates	PE-Cy7
Target Name	CD3E
Other Names	CD3E
Accession No.	Uniprot:P22646Gene ID:12501
Uniprot	P22646
GeneID	12501
Calculated MW	20
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only