

SDF1 beta Polyclonal Antibody FITC Conjugated

Catalog No: #C90313FITC



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	SDF1 beta Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-93/93
Conjugates	FITC
Target Name	SDF1 beta
Other Names	SDF1 beta
Accession No.	Uniprot:P48061Gene ID:6387
Uniprot	P48061
GeneID	6387
Calculated MW	8
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only