

TCP1 alpha Polyclonal Antibody Biotin Conjugated

Catalog No: #C90834Bio



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TCP1 alpha Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	225-340/556
Conjugates	Bio
Target Name	TCP1 alpha
Other Names	TCP1 alpha
Accession No.	Uniprot:P17987Gene ID:6950
Uniprot	P17987
GeneID	6950
Calculated MW	60
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only