



Caseine Kinase 1 alpha Polyclonal Antibody FITC Conjugated

Catalog No: #C90858FITC

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Caseine Kinase 1 alpha Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	221-337/337
Conjugates	FITC
Target Name	Caseine Kinase 1 alpha
Other Names	Caseine Kinase 1 alpha
Accession No.	Uniprot:P48729Gene ID:1452
Uniprot	P48729
GeneID	1452
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only