

Synaptopodin Polyclonal Antibody FITC Conjugated

Catalog No: #C91021FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Synaptopodin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	Synaptopodin
Other Names	Synaptopodin
Accession No.	Uniprot:Q8N3V7Gene ID:11346
Uniprot	Q8N3V7
GeneID	11346
Calculated MW	100
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only