

GM-CSFR alpha Polyclonal Antibody FITC Conjugated

Catalog No: #C91062FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GM-CSFR alpha Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-230/388
Conjugates	FITC
Target Name	GM-CSFR alpha
Other Names	GM-CSFR alpha
Accession No.	Uniprot:Q00941Gene ID:1438
Uniprot	Q00941
GeneID	1438
Calculated MW	38
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only