

SLC15A3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91554FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SLC15A3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	SLC15A3
Other Names	SLC15A3
Accession No.	Uniprot:Q8IY34Gene ID:51296
Uniprot	Q8IY34
GeneID	51296
Calculated MW	64
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt 1ug/Test

Note: This product is for in vitro research use only