

# Amphiregulin Polyclonal Antibody PE Conjugated

Catalog No: #C92065PE



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Amphiregulin Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	Amphiregulin
Other Names	Amphiregulin
Accession No.	Uniprot:P15514Gene ID:374
Uniprot	P15514
GeneID	374
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF 1:100-500

Note: This product is for in vitro research use only