

# GPR65 Polyclonal Antibody Bio Conjugated

Catalog No: #C92629Bio



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR65 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	GPR65
Other Names	GPR65
Accession No.	Uniprot:Q8IYL9Gene ID:8477
Uniprot	Q8IYL9
GeneID	8477
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000, IHC-P 1:100-500, IHC-F 1:100-500

Note: This product is for in vitro research use only