

IRE1 Polyclonal Antibody PE Conjugated

Catalog No: #C92713PE



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IRE1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	IRE1
Other Names	IRE1
Accession No.	Uniprot:O75460Gene ID:2081
Uniprot	O75460
GeneID	2081
Calculated MW	105
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only