

phospho-IRF3 (Ser386) Polyclonal Antibody PE Conjugated

Catalog No: #C92755PE

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	phospho-IRF3 (Ser386) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	S(p-S)LE
Conjugates	PE
Target Name	phospho-IRF3 (Ser386)
Other Names	phospho-IRF3 (Ser386)
Accession No.	Uniprot:Q14653Gene ID:3661
Uniprot	Q14653
GeneID	3661
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only