

# GAPDH-Loading Control mouse Monoclonal Antibody Biotin Conjugated

Catalog No: #C92848Bio



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GAPDH-Loading Control mouse Monoclonal Antibody Biotin Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	WB
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	GAPDH-Loading Control
Other Names	GAPDH-Loading Control
Accession No.	Uniprot:P04406Gene ID:2597
Uniprot	P04406
GeneID	2597
Calculated MW	38
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only