

EGFR mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92855FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EGFR mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	EGFR
Other Names	EGFR
Accession No.	Uniprot:P00533Gene ID:1956
Uniprot	P00533
GeneID	1956
Calculated MW	130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:100-500

Note: This product is for in vitro research use only