

# CD4 mouse Monoclonal Antibody APC Conjugated

Catalog No: #C92866APC



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD4 mouse Monoclonal Antibody APC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	IF
Species Reactivity	Hu Ms Rat
Conjugates	APC
Target Name	CD4
Other Names	CD4
Accession No.	Uniprot:P01730Gene ID:920
Uniprot	P01730
GeneID	920
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only