



alpha smooth muscle Actin mouse Monoclonal Antibody HRP Conjugated

Catalog No: #C92884HRP

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

| Description | |
|--------------------|---|
| Product Name | alpha smooth muscle Actin mouse Monoclonal Antibody HRP Conjugated |
| Host Species | mouse |
| Clonality | Monoclonal |
| Purification | affinity purified by Protein G |
| Applications | WB IHC ICC |
| Species Reactivity | Hu Ms Rat |
| Conjugates | HRP |
| Target Name | alpha smooth muscle Actin |
| Other Names | alpha smooth muscle Actin |
| Accession No. | Uniprot:P62736Gene ID:59 |
| Uniprot | P62736 |
| GeneID | 59 |
| Calculated MW | 42 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

| Application Details |
|---|
| WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 |

Note: This product is for in vitro research use only