



Cyclin B1 mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92899FITC

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cyclin B1 mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein G
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	Cyclin B1
Other Names	Cyclin B1
Accession No.	Uniprot:P14635Gene ID:891
Uniprot	P14635
GeneID	891
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only