



TGF beta 1 mouse Monoclonal Antibody FITC Conjugated

Catalog No: #C92902FITC

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TGF beta 1 mouse Monoclonal Antibody FITC Conjugated
Host Species	mouse
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	Hu Mo Rat
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	TGF beta 1
Other Names	TGF beta 1
Accession No.	Uniprot:P01137Gene ID:7040
Uniprot	P01137
GeneID	7040
Calculated MW	12.8/44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC IF

Note: This product is for in vitro research use only