

CD35 Rabbit mAb

Catalog No: #52204



Package Size: #52204-1 50ul #52204-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CD35 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S06-6D1
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human CD35
Conjugates	Unconjugated
Modification	Unmodification
Other Names	KN; C3BR; C4BR; CD35
Accession No.	Swiss-Prot:P17927GenelD:1378
Uniprot	P17927
GenelD	1378
Calculated MW	Calculated MW: 224 kDa; Observed MW: 280 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of CD35 in human tonsil lysates using CD35 antibody. Predicted band size: 224kDa. Observed band size: 280kDa.

Background

Swiss-Prot Acc.P17927.Mediates cellular binding of particles and immune complexes that have activated complement.(Microbial infection) Acts as a receptor for Epstein-Barr virus.

Note: This product is for in vitro research use only