

DHRS1 Rabbit mAb

Catalog No: #52252

Package Size: #52252-1 50ul #52252-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DHRS1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-9E4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human DHRS1
Conjugates	Unconjugated
Modification	Unmodification
Other Names	dehydrogenase/reductase 1;SDR19C1
Accession No.	Swiss-Prot:Q96LJ7GeneID:115817
Uniprot	Q96LJ7
GeneID	115817
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of DHRS1 in C6 cell lysates using DHRS1 Rabbit mAb(1:1000 diluted). Predicted band size:34kDa.Observed band size:34kDa.

Background

Swiss-Prot Acc.Q96LJ7.SitesFeature keyPosition(s)Description

Note: This product is for in vitro research use only