

Follistatin Rabbit mAb

Catalog No: #52287

Package Size: #52287-1 50ul #52287-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Follistatin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-8F4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human Follistatin
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene;
Accession No.	Swiss-Prot:P19883GeneID:10468
Uniprot	P19883
GeneID	10468
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of Follistatin in Rat Brain lysates using Follistatin Rabbit mAb(1:1000 diluted). Predicted band size:38kDa. Observed band size:38kDa.

Background

Swiss-Prot Acc.P19883.Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Note: This product is for in vitro research use only