

GAP43 Rabbit mAb

Catalog No: #52300



Package Size: #52300-1 50ul #52300-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GAP43 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S01-4F6
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human GAP43
Conjugates	Unconjugated
Modification	Unmodification
Other Names	B-50; PP46
Accession No.	Swiss-Prot:P17677GenelD:2596
Uniprot	P17677
GenelD	2596
Calculated MW	Calculated MW: 25 kDa; Observed MW: 46 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/1000; IHC: 1/50;

Images



Western blot detection of GAP43 in Rat Brain lysates using GAP43 Rabbit mAb(1:1000 diluted).Predicted band size:25kDa.Observed band size:46kDa.



Immunohistochemistry of GAP43 in paraffin-embedded Human Brain using GAP43 Rabbit mAb at dilution 1/50

Background

Swiss-Prot Acc.P17677.This protein is associated with nerve growth. It is a major component of the motile

Note: This product is for in vitro research use only